




Product Datasheet

BNIP3 antibody (orb31276)

Description	Rabbit polyclonal antibody to BNIP3
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IF, IHC-P
Immunogen	Synthetic Peptide
Target	BNIP3
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	21541
Uniprot ID	Q12983
NCBI	NP_004043.2
Dilution Range	IF/ICC: 1:50-100, WB: 1:1000, IHC-P: 1:50-100
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...



Confocal immunofluorescent analysis of m...