



Product Datasheet

PLA2R antibody (orb312775)

Description	PLA2R antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human PLA2R (1021-1100/1463 aa)
Target	PLA2R
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	153 kDa
Uniprot ID	Q13018
Dilution Range	WB=1:500-2000, IHC-P=1:100-500
Expiration Date	12 months from date of receipt.



Sample: Pancreas
(Mouse) Lysate at
40 ug...



Paraformaldehyde-
fixed, paraffin
embedde...