

Product Datasheet

POGZ antibody (orb312822)

Description	POGZ antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human POGZ/ZNF635 (501-600/1410 aa)
Target	POGZ/ZNF635
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	155 kDa
Uniprot ID	Q7Z3K3
Expiration Date	12 months from date of receipt.