


Product Datasheet

RGS12 antibody (orb312921)

Description	RGS12 antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC-Fr, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human RGS12 (601-700/1447 aa)
Target	RGS12
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	156 kDa
Uniprot ID	O14924
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of Human esophagus tissue.