

Product Datasheet

RGS16 (phospho-Tyr168) antibody (orb312924)

Description	RGS16 (phospho-Tyr168) antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr, IHC-P
Immunogen	KLH conjugated synthesised phosphopeptide derived from human RGS16 around the phosphorylation site of Tyr168 (DS(p-Y)PR)
Target	phospho-RGS16 (Tyr168)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	23 kDa
Uniprot ID	O15492
Expiration Date	12 months from date of receipt.