


Product Datasheet

RUNX2 antibody (orb312958)

Description	RUNX2 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	KLH conjugated synthetic peptide derived from human RUNX2 (201-300/521 aa)
Target	RUNX2
Preservatives	Preservative: 15mM Sodium Azide, Constituents: 1% rAlbumin, 0.01M PBS, pH 7.4
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	57(Hu)/67(Mo, Ra) kDa
Uniprot ID	Q13950
Dilution Range	Flow-Cyt=1µg/Test
Expiration Date	12 months from date of receipt.



Flow
cytometric
analysis of
U251 cell
us...