



Product Datasheet

Syncytin 1 antibody (orb312961)

Description	Syncytin 1 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from human Syncytin 1 (451-538/538 aa)
Target	Syncytin 1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	33/58 kDa
Uniprot ID	Q9UQF0
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.



Western blot:Lovo(Human)
Cell;Primary:An...



Western blot:Brain
(Mouse);Primary:Anti-
...