

Product Datasheet

SEC16A antibody (orb312992)

Description	This gene encodes a protein that forms part of the Sec16 complex. This protein has a role in protein...
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ICC, IF
Immunogen	KLH conjugated synthetic peptide derived from human SEC16A (1101-1200/2179 aa)
Target	SEC16A
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purified by Protein A
MW	Predicted Molecular Weight: 233kDa
Uniprot ID	O15027
Dilution Range	ICC=1:50-200, IF=1:50-200 - Optimal working dilutions must be determined by the end user.
Expiration Date	12 months from date of receipt.