


Product Datasheet

ZDHHC20 antibody (orb31317)

Description	ZDHHC20 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF
Immunogen	Recombinant Human Probable palmitoyltransferase ZDHHC20 protein (75-130AA)
Target	ZDHHC20
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q5W0Z9
Dilution Range	IF:1:50-200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of A549 cell...