

Product Datasheet

Mast Cell Tryptase antibody (orb313250)

Description	Mast Cell Tryptase antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-Fr
Immunogen	KLH conjugated synthetic peptide derived from human TPSD1 (171-242/242 aa)
Target	Mast Cell Tryptase
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	23 kDa
Uniprot ID	Q02844
Expiration Date	12 months from date of receipt.