

Product Datasheet

ERVK-6 antibody (FITC) (orb31341)

Description	ERVK-6 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Endogenous retrovirus group K member 6 Env polypeptide protein (90-632AA)
Target	ERVK-6
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q69384
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.