



Product Datasheet

Ubiquitin antibody (orb313417)


Description	Ubiquitin antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, ICC, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human Ubiquitin (63-74/76 aa)
Target	Ubiquitin
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	8.5 kDa
Uniprot ID	P0CG48
Dilution Range	WB=1:500-2000, IHC-P=1:100-500, ICC=1:100, Flow-Cyt=1ug/Test
Expiration Date	12 months from date of receipt.



Flow cytometry of Molt4.Primary Antibody...



ICC analysis of HepG2 cell using Ubiquit...



IHC-P analysis of human gastric carcinom...