



Product Datasheet

Vimentin (phospho-Ser55) antibody (orb313425)


Description	Vimentin (phospho-Ser55) antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	KLH conjugated synthesised phosphopeptide derived from human Vimentin around the phosphorylation site of Ser55 (YA(p-S)SP)
Target	phospho-Vimentin (Ser55)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	57 kDa
CAS Number	Vimentin
Uniprot ID	P08670
Dilution Range	WB=1:500-2000, Flow-Cyt=2ug/Test
Expiration Date	12 months from date of receipt.



Flow
cytometric
analysis of
A549 cells
I...



Western
blot analysis
of human
HT1080
ce...



Western
blot analysis
of human
hela cell...