

Product Datasheet

Streptococcus suis Suilysin antibody (orb313856)

Description	Streptococcus suis Suilysin antibody
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Full length SLY protein
Target	Streptococcus suis Suilysin
Preservatives	Preservative: 15mM Sodium Azide, Constituents: 1% rAlbumin, 0.01M PBS, pH 7.4
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Uniprot ID	T1R445
Expiration Date	12 months from date of receipt.