

Product Datasheet

CYP11A1 antibody (orb313866)

Description

CYP11A1 antibody


Species/Host

Rabbit

Reactivity


Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, IHC-P

Applications


IHC-P staining of
mouse ovarian tissue
u...

Immunogen

KLH conjugated synthetic peptide derived from human CYP11A1. Please contact us for the exact immunogen sequence. The peptide is available as orb374896.

Target

CYP11A1

Form/Appearance


10 mM PBS, 0.02% sodium azide

Concentration

- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl):
0.5 mg/ml

Storage

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.



Immunohistochemical
staining of mouse
ov...

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Purity

Polyclonal antibodies are purified by peptide affinity chromatography

MW

53 kDa

Uniprot ID

P05108, Q9QZ82, P14137

NCBI

000772, 000781

Entrez

1583

Dilution Range

IHC-P: 1:100-800

Expiration Date

12 months from date of receipt.