



Product Datasheet

ZO-1 antibody (orb313868)


Description	ZO-1 antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human ZO-1. Please contact us for the exact immunogen sequence. The peptide is available as orb374960.
Target	ZO-1
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
MW	191 kDa
Uniprot ID	P39447 , Q07157 , Q8R4T5
NCBI	003248 , 003257
Entrez	7082
Dilution Range	IF/ICC: 1:100-500, IHC-P: 1:100-500
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of mouse lung tissue using ZO-1 antibody.



IHC-P image of rat muscle tissue using ZO-1 antibody.



IHC-P staining of rat lymph node tissue using ZO-1 antibody.