



Product Datasheet

ZO-1 antibody (FITC) (orb313873)


Description	ZO-1 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	FITC
Tested Applications	ICC, IF
Immunogen	Synthetic peptide (please contact us for details)
Target	ZO-1 (FITC)
Form/Appearance	Liquid: PBS, 2 mM Sodium Azide, 30% glycerol
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
MW	191 kDa
Uniprot ID	P39447 , Q07157 , Q8R4T5
NCBI	NP_003248.3 , NM_003257.4
Entrez	7082
Dilution Range	IF/ICC:1:50-500 (based on 0.5 mg/ml)
Expiration Date	12 months from date of receipt.



IF analysis of rat colon tissue using an...



IF image of rat spleen tissue using anti...



Immunofluorescence analysis of rat splee...