

Product Datasheet

FAM182A antibody (FITC) (orb31434)

Description	FAM182A antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Protein FAM182A protein (1-154AA)
Target	FAM182A
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q5T1J6
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.