



Product Datasheet

CD3 epsilon antibody (orb314892)

Description	Hamster monoclonal antibody to CD3E
Species/Host	Hamster
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	FC, IP, WB
Immunogen	H-2Kb-specific murine cytotoxic T-lymphocyte (CTL) clone.
Target	CD3 epsilon
Preservatives	PBS with 0.02% Proclin 300.
Concentration	1 mg/ml
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Uniprot ID	P22646
Expiration Date	12 months from date of receipt.



WB analysis of Jurkat and Molt4 cell lines.



Flow-cytometric analysis of the dose-dependent binding of CD3 epsilon antibody to Jurkat cells.