



Product Datasheet

CD3 epsilon antibody (orb314895)

Description	Mouse monoclonal antibody to CD3E
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC
Immunogen	Human infant thymocytes.
Target	CD3 epsilon
Preservatives	PBS with 0.02% Proclin 300.
Concentration	1 mg/ml
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Application notes	This antibody binds to murine CD3 epsilon. It partially cross-blocks binding of clone 500A2 (Ab00212) indicating that overlapping epitopes are detected by these antibodies.
Isotype	IgG1
Clonality	Monoclonal
Purity	Purified
Uniprot ID	P07766
Expiration Date	12 months from date of receipt.



Flow-cytometric analysis of the dose-dep...



WB analysis of Jurkat and Molt4 cell lin...