


Product Datasheet

Human ADAM9 ELISA Kit (orb314899)

Description	Sandwich ELISA for quantitative detection of Human ADAM9 concentrations.
Reactivity	Human
Range	15.6 pg/ml - 1000 pg/ml
Tested Applications	ELISA
Target	ADAM9
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Experiment Notes: Cell Culture Supernatants, Serum, Plasma
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	Q13443
Sensitivity	10 pg/ml
Entrez	8754
Expiration Date	Please enquire.



Representative
Standard
Curve