


Product Datasheet

Human ADAMTS1 ELISA Kit (orb314900)

Description	Sandwich ELISA for quantitative detection of Human ADAMTS1 concentrations.
Reactivity	Human
Range	93.7 pg/ml - 6000 pg/ml
Tested Applications	ELISA
Target	ADAMTS1
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Experiment Notes: Cell Culture Supernatants, Serum, Plasma
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	Q9UHI8
Sensitivity	10 pg/ml
Entrez	9510
Expiration Date	Please enquire.



Representative
Standard
Curve