

Product Datasheet

Human CD321 ELISA Kit (orb314912)

Description Sandwich ELISA for quantitative detection of Human CD321 concentrations.

Reactivity Human

Range 46.9 pg/ml - 3000 pg/ml

Tested Applications ELISA

Target F11R

Storage Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.

Note For research use only

Application notes Experiment Notes: Cell Culture Supernatants, Serum


Sample Types Cell Culture Supernates, Serum

Uniprot ID [Q9Y624](#)

Sensitivity 10 pg/ml

Entrez [50848](#)

Expiration Date Please enquire.



Representative
Standard
Curve