


Product Datasheet

Human CD32 ELISA Kit (orb314913)

Description	Sandwich ELISA for quantitative detection of Human CD32 concentrations.
Reactivity	Human
Range	46.9 pg/ml - 3000 pg/ml
Tested Applications	ELISA
Target	FCGR2B; FCGR2C
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Experiment Notes: Cell Culture Supernatants, Serum, Plasma
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	P31994, P31995
Sensitivity	1 pg/ml
Entrez	9103, 2213
Expiration Date	Please enquire.



Representative
Standard
Curve