


Product Datasheet

Human IL-37 ELISA Kit (orb314920)

Description	Sandwich ELISA for quantitative detection of Human IL-37 concentrations.
Reactivity	Human
Range	31.2 pg/ml - 2000 pg/ml
Tested Applications	ELISA
Target	IL37
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Experiment Notes: Cell Culture Supernatants, Serum, Plasma
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	Q9NZH6
Sensitivity	15 pg/ml
Entrez	27178
Expiration Date	Please enquire.



Representative
Standard
Curve