


Product Datasheet

Human Mesothelin ELISA Kit (orb314923)

Description	Human Mesothelin ELISA Kit
Reactivity	Human
Range	156 pg/ml - 10000 pg/ml
Tested Applications	ELISA
Target	MSLN
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Experiment Notes: Cell Culture Supernatants, Serum, Plasma, Urine
Sample Types	Cell Culture Supernates, Serum, Plasma, Urine
Uniprot ID	Q13421
Sensitivity	10 pg/ml
Entrez	10232
Expiration Date	Please enquire.



Representative
Standard
Curve