


Product Datasheet

Human CD304 ELISA Kit (orb314925)

Description	Sandwich ELISA for quantitative detection of Human CD304 concentrations.
Reactivity	Human
Range	780 pg/ml - 50000 pg/ml
Tested Applications	ELISA
Target	NRP1
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Experiment Notes: Cell Culture Supernatants, Serum, Plasma
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	O14786
Sensitivity	20 pg/ml
Entrez	8829
Expiration Date	Please enquire.



Representative
Standard
Curve