


Product Datasheet

Human Syndecan 1 ELISA Kit (orb314927)

Description	Human Syndecan 1 ELISA Kit
Reactivity	Human
Range	156 pg/ml - 10000 pg/ml
Tested Applications	ELISA
Target	SDC1
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Experiment Notes: Cell Culture Supernatants, Serum
Sample Types	Cell Culture Supernates, Serum, Plasma, Cell Lysates, Tissue Homogenates
Uniprot ID	P18827
Sensitivity	10 pg/ml
Entrez	6382
Expiration Date	Please enquire.



Representative
Standard
Curve