

Product Datasheet

Human Hevin ELISA Kit (orb314932)

Description Sandwich ELISA for quantitative detection of Human Hevin concentrations.

Reactivity Human

Range 46.9 pg/ml - 3000 pg/ml

Tested Applications ELISA

Target SPARCL1

Storage Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.

Note For research use only

Application notes Experiment Notes: Cell Culture Supernatants, Serum


Sample Types Cell Culture Supernates, Serum

Uniprot ID **Q14515**

Sensitivity 10 pg/ml

Entrez **8404**

Expiration Date Please enquire.



Representative
Standard
Curve