


Product Datasheet

Mouse Angiopoietin-1 ELISA Kit (orb314936)

Description	Sandwich ELISA for quantitative detection of Mouse Angiopoietin-1 concentrations.
Reactivity	Mouse
Range	156 pg/ml - 10000 pg/ml
Tested Applications	ELISA
Target	ANGPT1
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Experiment Notes: Cell Culture Supernatants, Serum, Plasma
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	008538
Sensitivity	10 pg/ml
Entrez	11600
Expiration Date	Please enquire.



Representative
Standard
Curve