


Product Datasheet

Mouse CD144 Cadherin ELISA Kit (orb314944)

Description	Sandwich ELISA for quantitative detection of Mouse CD144 Cadherin concentrations.
Reactivity	Mouse
Range	312 pg/ml - 20000 pg/ml
Tested Applications	ELISA
Target	CDH5
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Experiment Notes: Cell Culture Supernatants, Serum, Plasma
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	P55284
Sensitivity	10 pg/ml
Entrez	12562
Expiration Date	Please enquire.



Representative
Standard
Curve