


Product Datasheet

Mouse FGF21 ELISA Kit (orb314948)

Description	Sandwich ELISA for quantitative detection of Mouse FGF21 concentrations.
Reactivity	Mouse
Range	15.6 pg/ml - 1000 pg/ml
Tested Applications	ELISA
Target	FGF21
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Experiment Notes: Cell Culture Supernatants, Serum, Plasma
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	Q9JJN1
Sensitivity	8 pg/ml
Entrez	56636
Expiration Date	Please enquire.



Representative
Standard
Curve