


Product Datasheet

Mouse IL-17D ELISA Kit (orb314953)

Description	Sandwich ELISA for quantitative detection of Mouse IL-17D concentrations.
Reactivity	Mouse
Range	62.5 pg/ml - 4000 pg/ml
Tested Applications	ELISA
Target	IL17D
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Experiment Notes: Cell Culture Supernatants, Serum, Plasma
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	Q8K4C4
Sensitivity	10 pg/ml
Entrez	239114
Expiration Date	Please enquire.



Representative
Standard
Curve