

Product Datasheet

Mouse CD127 ELISA Kit (orb314956)

Description Sandwich ELISA for quantitative detection of Mouse CD127 concentrations.

Reactivity Mouse

Range 31.2 pg/ml - 2000 pg/ml

Tested Applications ELISA

Target IL7R

Storage Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.

Note For research use only

Application notes Experiment Notes: Cell Culture Supernatants, Serum


Sample Types Cell Culture Supernates, Serum

Uniprot ID **P16872**

Sensitivity 10 pg/ml

Entrez **16197**

Expiration Date Please enquire.



Representative
Standard
Curve