


Product Datasheet

Mouse PLTP ELISA Kit (orb314961)

Description	Sandwich ELISA for quantitative detection of Mouse PLTP concentrations.
Reactivity	Mouse
Range	780 pg/ml - 50000 pg/ml
Tested Applications	ELISA
Target	PLTP
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Experiment Notes: Cell Culture Supernatants, Serum, Plasma
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	P55065
Sensitivity	20 pg/ml
Entrez	18830
Expiration Date	Please enquire.



Representative
Standard
Curve