


Product Datasheet

Rat Adiponectin ELISA Kit (orb314971)

Description	Sandwich ELISA for quantitative detection of Rat Adiponectin concentrations.
Reactivity	Rat
Range	156 pg/ml - 10000 pg/ml
Tested Applications	ELISA
Target	ADIPOQ
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Experiment Notes: Cell Culture Supernatants, Serum, Plasma, Tissue Homogenates, Urine
Sample Types	Cell Culture Supernates, Serum, Plasma, Tissue Homogenates, Urine
Uniprot ID	Q8K3R4
Sensitivity	10 pg/ml
Entrez	246253
Expiration Date	Please enquire.



Representative
Standard
Curve