


Product Datasheet

Rat CD325 ELISA Kit (orb314976)

Description	Sandwich ELISA for quantitative detection of Rat CD325 concentrations.
Reactivity	Rat
Range	780 pg/ml - 50000 pg/ml
Tested Applications	ELISA
Target	CDH2
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Experiment Notes: Cell Culture Supernatants, Serum, Plasma
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	Q9Z1Y3
Sensitivity	20 pg/ml
Entrez	83501
Expiration Date	Please enquire.



Representative
Standard
Curve