


Product Datasheet

Rat VEGFR3 ELISA Kit (orb314979)

Description	Sandwich ELISA for quantitative detection of Rat VEGFR3 concentrations.
Reactivity	Rat
Range	156 pg/ml - 10000 pg/ml
Tested Applications	ELISA
Target	FLT4
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Experiment Notes: Cell Culture Supernatants, Serum
Sample Types	Cell Culture Supernates, Serum
Uniprot ID	Q91ZT1
Sensitivity	10 pg/ml
Entrez	114110
Expiration Date	Please enquire.



Representative
Standard
Curve