


Product Datasheet

Human IL-13 ELISA Kit (orb314988)

Description	Sandwich ELISA for quantitative detection of Human IL-13 concentrations.
Reactivity	Human
Range	93.8 pg/ml - 6000 pg/ml
Tested Applications	ELISA
Target	IL13
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Experiment Notes: Cell Culture Supernatants, Serum, Plasma, Tissue Homogenates
Sample Types	Cell Culture Supernates, Serum, Plasma, Tissue Homogenates
Uniprot ID	P35225
Sensitivity	40 pg/ml
Entrez	3596
Expiration Date	Please enquire.



Representative
Standard
Curve