


Product Datasheet

Human IFN gamma ELISA Kit (orb314998)

Description	Human IFN gamma ELISA Kit
Reactivity	Human
Range	15.6 pg/ml - 1000 pg/ml
Tested Applications	ELISA
Target	IFNG
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Experiment Notes: Cell Culture Supernatants, Serum, Plasma, Tissue Homogenates, Urine
Sample Types	Cell Culture Supernates, Serum, Plasma, Tissue Homogenates, Urine
Uniprot ID	P01579
Sensitivity	4 pg/ml
Entrez	3458
Expiration Date	Please enquire.



Representative
Standard
Curve