


Product Datasheet

Human TNF alpha ELISA Kit (orb314999)

Description	Sandwich ELISA for quantitative detection of Human TNF alpha concentrations.
Reactivity	Human
Range	15.6 pg/ml - 1000 pg/ml
Tested Applications	ELISA
Target	TNF
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Experiment Notes: Cell Culture Supernatants, Serum, Plasma, Tissue Homogenates, Urine
Sample Types	Cell Culture Supernates, Serum, Plasma, Tissue Homogenates, Urine
Uniprot ID	P01375
Sensitivity	7 pg/ml
Entrez	7124
Expiration Date	Please enquire.



Representative
Standard
Curve