


Product Datasheet

Human IL-32 ELISA Kit (orb315018)

Description	Sandwich ELISA for quantitative detection of Human IL-32 concentrations.
Reactivity	Human
Range	7.8 pg/ml - 500 pg/ml
Tested Applications	ELISA
Target	IL32
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Experiment Notes: Cell Culture Supernatants, Serum, Plasma
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	P24001
Sensitivity	4 pg/ml
Entrez	9235
Expiration Date	Please enquire.



Representative
Standard
Curve