

Product Datasheet

Human IL-31RB ELISA Kit (orb315024)

Description Sandwich ELISA for quantitative detection of Human IL-31RB concentrations.

Reactivity Human

Range 312 pg/ml - 20000 pg/ml

Tested Applications ELISA

Target OSMR

Storage Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.

Note For research use only

Application notes Experiment Notes: Cell Culture Supernatants, Serum, Plasma


Sample Types Cell Culture Supernates, Serum, Plasma

Uniprot ID **Q99650**

Sensitivity 80 pg/ml

Entrez **9180**

Expiration Date Please enquire.



Representative
Standard
Curve