


Product Datasheet

Human Glucagon ELISA Kit (orb315026)

Description	Human Glucagon ELISA Kit
Reactivity	Human
Range	62.5 pg/ml - 4000 pg/ml
Tested Applications	ELISA
Target	GCG
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Experiment Notes: Cell Culture Supernatants, Serum, Plasma
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	P01275
Sensitivity	30 pg/ml
Entrez	2641
Expiration Date	Please enquire.



Representative
Standard
Curve