

## Product Datasheet

### Human GDF15 ELISA Kit (orb315035)

**Description** Sandwich ELISA for quantitative detection of Human GDF15 concentrations.

**Reactivity** Human

**Range** 7.8 pg/ml - 500 pg/ml

**Tested Applications** ELISA

**Target** GDF15

**Storage** Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.

**Note** For research use only

**Application notes** Experiment Notes: Cell Culture Supernatants, Serum, Plasma


**Sample Types** Cell Culture Supernates, Serum, Plasma

**Uniprot ID** **Q99988**

**Sensitivity** 4 pg/ml

**Entrez** **9518**

**Expiration Date** Please enquire.



Representative  
Standard  
Curve