


Product Datasheet

Mouse IFN gamma ELISA Kit (orb315037)

Description	Mouse IFN gamma ELISA Kit
Reactivity	Mouse
Range	7.8 pg/ml - 500 pg/ml
Tested Applications	ELISA
Target	IFNG
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Experiment Notes: Cell Culture Supernatants, Serum, Plasma, Tissue Homogenates
Sample Types	Cell Culture Supernates, Serum, Plasma, Tissue Homogenates
Uniprot ID	P01580
Sensitivity	4 pg/ml
Entrez	15978
Expiration Date	Please enquire.



Representative
Standard
Curve