


Product Datasheet

Mouse IL-5 ELISA Kit (orb315047)

Description	Sandwich ELISA for quantitative detection of Mouse IL-5 concentrations.
Reactivity	Mouse
Range	15.6 pg/ml - 1000 pg/ml
Tested Applications	ELISA
Target	IL5
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Experiment Notes: Cell Culture Supernatants, Serum, Plasma, Tissue Homogenates
Sample Types	Cell Culture Supernates, Serum, Plasma, Tissue Homogenates
Uniprot ID	P04401
Sensitivity	2 pg/ml
Entrez	16191
Expiration Date	Please enquire.



Representative
Standard
Curve