


Product Datasheet

Mouse IL12 p40 ELISA Kit (orb315048)

Description	Sandwich ELISA for quantitative detection of Mouse IL12 p40 concentrations.
Reactivity	Mouse
Range	7.8 pg/ml - 500 pg/ml
Tested Applications	ELISA
Target	IL12B
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Experiment Notes: Cell Culture Supernatants, Serum, Plasma, Tissue Homogenates
Sample Types	Cell Culture Supernates, Serum, Plasma, Tissue Homogenates
Uniprot ID	P43432
Sensitivity	4 pg/ml
Entrez	16160
Expiration Date	Please enquire.



Representative
Standard
Curve