


Product Datasheet

Mouse GM-CSF ELISA Kit (orb315050)

Description	Sandwich ELISA for quantitative detection of Mouse GM-CSF concentrations.
Reactivity	Mouse
Range	31.2 pg/ml - 2000 pg/ml
Tested Applications	ELISA
Target	CSF2
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Experiment Notes: Cell Culture Supernatants, Serum, Plasma
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	P01587
Sensitivity	15 pg/ml
Entrez	12981
Expiration Date	Please enquire.



Representative
Standard
Curve