


Product Datasheet

Mouse TSLP ELISA Kit (orb315051)

Description	Sandwich ELISA for quantitative detection of Mouse TSLP concentrations.
Reactivity	Mouse
Range	15.6 pg/ml - 1000 pg/ml
Tested Applications	ELISA
Target	TSLP
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Experiment Notes: Cell Culture Supernatants, Serum, Plasma
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	Q9JIE6
Sensitivity	7 pg/ml
Entrez	53603
Expiration Date	Please enquire.



Representative
Standard
Curve